PROGRAMIX 25-50-100

UGIN'
dentaire

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INSTRUCTIONS FOR USE

INTRODUCTION	18
SPECIFICATION	19
INSTALLATION AND SET -UP	20
DESCRIPTION OF CONTROL STATION	21
PROGRAMMING /CALIBRATION	22
PROTECTIONS	27
TROUBLE SHOOTING	28
REPLACING THE PYROMETER	29
REPLACING THE MUFLE	29
°C/°F CONVERSION CHART	30
TABLE OF SINGLE PROGRAMS	31
TABLE OF DOUBLE PROGRAMS	32

Translation into English from original in French

INTRODUCTION

Technological advances in investment materials have put increasing demands on the performance of the furnaces used.

Accurate casting of a given investment material no longer depends on the precision of the firing cycle alone.

The PROGRAMIX furnace, manufactured using the latest heating and control techniques, will offer you invaluable assistance in achieving top-quality results.


We recommend that you read this instruction manual attentively before beginning to operate your furnace.

SPECIFICATION

	PX 25	PX 50	PX 100
- Maximum temperature - Continuous operation - Peak operation	1050°C/ 1150°C/	1050°C/ 1150°C/	1050°C/ 1150°C/
Timer programming	0 to 100 h	0 to 100 h	0 to 100 h
Regulated heat rise rate	1 to 15°C	1 to 15°C	1 to 15°C
Heat rise rate Full position	16°C and +	16°C and +	16°C and +
Number of programs	30	30	30
Program settings 0 through 19	3 stages	3 stages	3 stages
Program settings 20 through 29	6 stages	6 stages	6 stages
Power requirement	230 V 50-60	230V 50-60	230 V 50-60
Rated power	1400 W	2000 W	3000 W
Connection of exhauster flue	230V	230 V	230V
- Furnace dimensions - Height - Width - Depth	460 mm 350 mm 370 mm	500 mm 390 mm 440 mm	590 mm 465 mm 555 mm
- Muffle dimensions - Height - Width - Depth	100 mm 150 mm 100 mm	120 mm 170 mm 210 mm	150 mm 220 mm 300 mm
Number of casting rings (dia. 80.length 55)	2 à 3	4 à 6	8 à 12
Weight	28 kg/62Lbs	34 kg/75Lbs	56 kg/124Lbs

INSTALLATION AND SET- UP

1. Vent outlet
2. Fuse holder
3. Power socket
4. Flue exhauster socket
5. Identification plate
6. Location of TRIAC
7. Thermocouple cover

► After unpacking the machine, verify its condition and notify the shipper and supplier of any damage.

► Remember to have the deliveryman sign the shipper order on which any reservations have been noted and to confirm the latter to the shipper by registered letter with acknowledgement of receipt.

INSTALLATION FO THE FLUE EXHAUSTER: (available as an option)

The flue exhauster is the normal complement of the pre-drying furnace.
It is connected electrically to the socket **4** and mechanically in place of the vent outlet **1**.

The flue exhauster handles the obvious job of removing waste gases from the laboratory.
It also helps keep the front panel and controls clean at an times.

INSTALLATION OF THE FURNACE:


► After checking that the mains voltage is the same as that indicated on the identification plate **5** connect the cord to the socket **3**.

► Place the furnace in a ventilated area at a distance of at least 10 cm from any other vertical surfaces.

► Then connect the power cord to a standard 10/16A socket with a grounding terminal.


► Flip the master switch located on the right and side of the furnace.

DESCRIPTION OF THE CONTROL STATION

- 1.** Panel of operating parameters.
- 2.** Program number display (0 to 29).
- 3.** Temperature/time display.
- 4.** Indicator of the unit displayed on **3**.


- 5.** Mode key
- 6.** Forward key
- 7.** Forward key
- 8.** STOP or START touch control
- 9.** "CYCLE" advance key
- 10.** Program selection key
- 11.** Indicators:


- Double (used for programs 20 through 29)

- Calibration

- Heating element

- Ventilation

PROGRAMMING

1) PANEL OF OPERATING PARAMETERS

- To program the operating parameters, press key .
- When the program number on display **2** begins blinking, you may enter the data by pressing the keys  .

To switch to another parameter, use the  key.


When you have programmed all the parameters, "lock" the program by pressing key  again.

Explanation of the parameters:

1.1) Start time - delay. 0 to 99h59min.

This timer function is used to delay the beginning of the firing cycle.

A Count-down before start time is displayed on **3**.

Should you not wish to use this function, set the parameter at 0h 0 min and heating will begin immediately (or use the  key after pressing the  touch control)

1^e STAGE

1.2) 1st heating rate (Rampe N° 1)

The heating rate is the heat rise per unit time (rate of rise).

There are three positions:

- **STOP position** - the furnace stops heating.
This is only used when the next stage is at a lower temperature.
- **Position 1 to 15°/mn**
Heating is regulated and will not rise above the programmed value. However, it may be lower than that value, especially in the higher temperature ranges or if the muffle is well loaded.
- **"FULL" Position**
This position corresponds to the maximum furnace out put. Here, the heat rise rate is no longer regulated.

1.3) Température 0 to 1200°C

Temperature setting for the 1 st stage.

1.4) Hold time 0 to 240 mn

Hold time for the temperature of the 1 st stage.

2^e STAGE

1.5) Heating rate 2 (RP2) 1 to 15° /min - FULL - STOP

Heat rise to reach the second stage. The principle is the same as for heating rate 1, except that the programmed values may differ.

1.6) Temperature of the 2nd stage - 0 to 1200.C (see paragraph 3 above)

1.7) Hold time. 0 to 240 min (see paragraph 4 above)

3^e STAGE

1.8) Heating rate 3 (RP3) - 1 to 15°C/ min - FULL - STOP

1.9) Temperature 0 to 1200° C

2.0) Hold time. 0 to 240 min

2.1) Final time-delay 0 to 240 min

This timer function sets a time delay before the firing cycle is stopped completely.

Once you have programmed your full cycle, you may program the furnace to hold at the last temperature for 0 to 240 min, in addition to the hold time.

This means that you can choose the time you need to perform your castings.

2) PROGRAM NUMBER DISPLAY 0 TO 29 :

If the control station is on stand-by, the program number can be changed using the keys  

The control station has a storage capacity of 30 programs.

Programs 0 through 19 are 3-stage programs.

Programs 20 through 29 are 6-stage programs.

 The program number cannot be changed while the furnace is in operation (the program number blinks, key  or when the cycle has begun ( control).

3) TEMPERATURE /TIME DISPLAY

This display gives you the value of the parameters while setting up the program

During the heating cycle, you can choose to display either:

- The temperature during the heat rise and the countdown during the hold periods,
- Or the temperature only.

In that case, press the  key.

A dot then appears at the far right of the display **3** to indicate that you have chosen the temperature mode.

4) INDICATOR OF THE UNIT DISPLAYED ON 3

Two red lights indicate units of temperature and heating rate ($^{\circ}\text{C}$ and $^{\circ}\text{C} / \text{min}$) and two green lights indicate units of time (hours and minutes).

5) MODE KEY

This key is pressed to switch from the temperature/time display mode to the temperature only display mode. See paragraph 3 above.

6) FORWARD KEY

7) REVERSE KEY

7.1) In stand-by position:

Pressing the  or the  key modifies the program number.

7.2) In programming position:

Pressing the  key increases the value of the parameter displayed on **3**.

Pressing the key  decreases that value.

8) STOP OR START TOUCH CONTROL

After checking the number of the program selected, the  control is pressed to start the heating cycle. All the indicators on the panel light up simultaneously.

By pressing the  touch control a second time, the cycle is cancelled.

9) CYCLE ADVANCE KEY

During program set-up, this key is pressed to choose the parameters to be programmed.

.During the heating cycle, the  key is pressed to skip a stage and move on to the next one.

10) PROGRAM SELECTION KEY

This key accesses program selection.

Once  is pressed, program parameters can be entered or modified using the   or controls and the  key.

During that time, the program number blinks.

Once the program has been entered, pressing the  key a second time "locks" the program. The program number then stops blinking.

11) INDICATORS

► Double


- Programs 0 through 19 are simple three-stage programs.
- Programs 20 through 29 are double 6-stage programs.

While programs 20 through 29 are being set up, this indicator goes off for stages 1 to 3 and on for stages 4 to 6.

► Calibration


The "calibration" function is used to adjust the heating temperature in the furnace within a range of + or – 50°C, by increments of 10°C (see the following page)

► Heating element

This indicator visualizes the electric pulses transmitted to the muffle. It continues to blink during the heat rise or hold time.

► Ventilation

A flue exhauster (available as an option) may be connected via the socket located on the back panel of the PROGRAMIX furnaces (25-50-100).

It is programmed after the final time-delay parameter has been set by pressing the  key and then  .

The exhauster starts operating at the beginning of the heating cycle and stops when the programmed temperature is reached.

TRIAL RUN:

After installing your furnace according to the instructions, on page SPECIFICATION, flip the master switch located on the right side of the furnace.

The main display indicates the temperature at 3 and the program number at 2.

■ **To change the program number, use**  .

■ **To visualize or modify the data in a program:**

 >   >  >   > etc

(during this operation, the program number blinks).

■ **To "lock" the program:**

Press  second time.

(The program number stops blinking)

■ To Launch the cycle:

Press the  touch control, all the panel indicators light up.

Each of these indicators will go out at the appropriate time during the heating cycle.

N.B.:

 .To move the program manually from one stage to the next, use the  key.

 . To cancel an active cycle, press the  touch control.

CALIBRATION

All UGIN'DENTAIRE furnaces have been subjected to an IPTS 68 (International Practical Temperature Scale of 1968) temperature test using a temperature calibrator.

However, should you not find the initial setting satisfactory; you can test the temperature using a pure silver wire. Pure silver melts at 960°C.

TEST PROCEDURE

-Place the silver wire on a ceramic support in the middle of the muffle.

- Raise the temperature by stages. starting at 940°C with a 5 min hold, then 950°C/5 min, then 960°C/ 5 min and, If need be, 970°C/min, until the test silver melts. (Open the furnace door briefly after each stage to check on the silver wire).

Let us assume, for example, that the wire melted at 950°C. In that case, the calibration should be modified as follows:

Press  and then 

The figure **5** appears on the main display and the "Calib" indicator lights up.

Using the  key, raise the position number to **6**.

Press  again.

N.B:

 .Each position corresponds to a temperature change of 10° C relative to the starting point. Hence, it is normal to observe that same deviation from ambient temperature when the furnace is cold.

PROTECTIONS

PROGRAM PROTECTION:

- The control station is equipped with a lithium battery to keep the programs stored in the memory
- for approximately five years.

After that period, we recommend that the battery be replaced.

OVERVOLTAGE PROTECTION:

The power supply to a furnace microprocessor must be perfectly filtered.

The PROGRAMIX furnaces are protected against over voltage by a system of active and passive filters.

Nevertheless. Abnormally high over voltages may cause some interference in the programs. As a safety precaution, we advise you to record your programs on the table prepared for that purpose (see the following page)

POWER FLUCTUATION:

Fluctuations in the mains voltage (within a range of $\pm 20V$) have no impact on the operation of the temperature programmer. However, the heat rise may be lower than the programmed value if the line drop continues.


POWER FAILURE:

The active cycle is interrupted.

TROUBLE SHOOTING

The **PROGRAMIX 25- 50-100** furnace are remarkably reliable. However, despite all the care taken in producing it, a failure is always possible, as with any other piece of equipment.

IF THE FURNACE DOES NOT HEAT

- **The main display shows "dF sond"** (thermocouple failure)
Replace the pyrometer as described (see the following page)
- The  **indicator stays lit continuously, but there is no temperature rise.**
 - Remove the lower rear panel and check resistance at the muffle with a multi-meter.

PROGRAMIX 25	⇒	35 ohms
PROGRAMIX 50	⇒	24 ohms
PROGRAMIX 100	⇒	15 ohms

- **If the resistor is damaged**
 - Replace the muffle as described (see the following page).

THE FURNACE CONTINUES TO HEAT BEYOND PROGRAMMED LIMITS

- Replace the TRIAC located on the rear panel.

REPLACING THE PYROMETER

- ▶ With a screwdriver, unscrew the thermocouple cover located at the rear of the furnace near the vent.
- ▶ Then unscrew the pyrometer and disconnect the two wires attached to it.
- ▶ Screw in the new pyrometer and re-connect the two wires, making sure to hook up positive and negative poles correctly:
 - + **Green wire**
 - **White wire**

Attention: the negative strand on the pyrometer is slightly magnetic.

- ▶ Let the furnace heat at a steady 200°C to check that it is working properly.

REPLACING THE MUFFLE

- ▶ Flip the master switch to off. Disconnect the power cord. After dismantling the pyrometer, unscrew the rear panel.
- ▶ Remove the insulation plates. Disconnect the muffle. Take the fibre insulation off the upper half. Pull the muffle out from the rear.
- ▶ Install the new muffle, taking care to keep the resistor terminals at least 5 cm away from the heating element at all times.
- ▶ Reassemble the various components by reversing the dismantling order.
- ▶ Dehydrate the furnace by heating it at 200°C for one hour. It is then ready for use again.

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